BioGlue®

Instructions for Use

Surgical Adhesive

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CE0088

ENGLISH

Read *Instructions for Use* prior to using this product.

RxONLY = Prescription Device.

In the United States of America, Federal Law restricts a prescription device to sale by or on the order of a physician.

In all other regions of the world, this device should not be used except on the order of a physician.

NOT STERILE = Non-Sterile

The BioGlue Surgical Adhesive (BioGlue) Reusable Syringe Delivery Device is to be used in conjunction with 5mL or 10mL BioGlue syringes and applicator tips that are supplied separately. Refer to the BioGlue Surgical Adhesive *Instructions for Use* included with the BioGlue syringes for additional information such as indications for use, contraindications, warnings, and precautions.

HOW SUPPLIED

The BioGlue Reusable Syringe Delivery Device is supplied non-sterile and must be sterilized prior to its first use. The Delivery Device may then be resterilized for reuse nine additional times. Retain these *Instructions for Use* for future reference.

MAINTAINING DELIVERY DEVICE EFFECTIVENESS

Precautions

1. Do not use glutaraldehyde-based cleaners or other cleaners corrosive to metal parts as these may corrode parts of the Delivery Device and impede its function.
2. Do not use ultrasonic cleaning methods as these may lead to material fatigue and to failure of the Delivery Device.
3. Do not use dry heat or ethylene oxide sterilization methods.
4. Do not clean with enzyme-based detergents since these may reduce the number of times the Delivery Device may be reused.

Cleaning

Inspect the Delivery Device for signs of deterioration or damage such as cracking, blistering, or pitting. Discard if deterioration or damage is observed. Remove the plunger from the handle by lifting the plunger latch upward and drawing the plunger forward. The Delivery Device may be satisfactorily cleaned by thoroughly washing with a mild detergent according to the following directions:

Step	Time	Cleaning Criteria
1	As needed	Remove the plunger from the delivery device by lifting the plunger latch upward and drawing the plunger forward.
2	2 minutes	Rinse with tap water prior to soak.
3	5 minutes	Place the Delivery Device in a 10 - 40°C (50 - 104°F) water bath containing Alconox (diluted according to manufacturer's instructions). Clean all crevices using a soft bristle brush to remove any visible debris.
4	2 minutes	Rinse and flush with DI water.

Sterilization

Sterilize by steam autoclaving only. The Delivery Device should first be cleaned and placed in appropriate packing materials. The following autoclave settings have been tested to ensure sterility and product functionality:

Cycle Temperature (C)	Cycle Temperature (F)	Cycle Time	Cycle Type
121°C	250°F	90 minutes	Pre-vacuum
128-132°C	262-270°F	24 minutes	Pre-vacuum
134°C	273°F	3-24 minutes	Pre-vacuum
132°C	270°F	10 minutes	Gravity

Exposures outside of these settings have not been evaluated. The Delivery Device is supplied non-sterile and must be sterilized prior to its first use. The Delivery Device may then be cleaned and reesterilized for reuse nine additional times under normal use conditions.

The user must validate the equipment chosen to sterilize the device. The effectiveness of each sterilization cycle should be monitored using biological indicators and/or other means.

NOTE: The user should verify that the Delivery Device is still functional after each cleaning and sterilization. This is accomplished by assembling the Delivery Device and observing that the plunger moves forward with each depression of the lever.

DELIVERY DEVICE PREPARATION

The BioGlue Surgical Adhesive and its Reusable Syringe Delivery Device consists of: main delivery unit, black Delivery Device plunger, syringe, and applicator tips.

1. Remove the main delivery unit and black Delivery Device plunger from their packaging. While pushing upward on the small latch at the back of the main unit, insert the plunger, ribbed side downward, into the front of the main unit.

NOTE: The white plunger supplied with the syringe will not work with the Delivery Device.

2. Still holding the small latch upward, retract the Delivery Device plunger completely back.
3. Remove a syringe from its packaging. Align the large and small barrels of the syringe with the corresponding Delivery Device plunger heads in the main delivery unit. From the side of the Delivery Device, slide the syringe into position, leading with the large barrel. Slide the Delivery Device plunger forward until resistance is felt.

CAUTION: Do not compress the main delivery unit trigger mechanism while attaching the syringe.

4. Remove an applicator tip from its packaging and inspect the collar portion of the tip to ensure that the triangular notch is directly over the larger port. If not, rotate the locking collar on the shaft until the notch is over the larger port.
5. While firmly grasping the syringe, nose upward, turn the syringe cap 90° counterclockwise and remove the cap by rocking it from side-to-side. Align the tip with the syringe using the corresponding notches on each and place the tip on the syringe.

CAUTION: Take care not to spill solution from the syringe during assembly.

CAUTION: Do not compress the main delivery unit trigger mechanism while attaching the applicator tip.

ENGLISH

- Lock the applicator tip in place by pushing the tip firmly toward the syringe and rotating the tip collar 90° clockwise.

NOTE: Before using BioGlue in the procedure, the residual air space must be removed and the applicator tip primed. Refer to [Site Preparation, Syringe Air Space Removal and Applicator Tip Priming](#).

- If using an applicator tip with a flexible extension, a desired angle may be created by bending the extension at the appropriate location to the desired angle and holding for 3-5 seconds. The angle created should be maintained for up to 5 minutes.
- To remove occluded applicator tips, grasp the applicator tip collar, rotate the tip collar counterclockwise, and lift the tip off the syringe by rocking it side-to-side.

SITE PREPARATION, SYRINGE AIR SPACE REMOVAL AND APPLICATOR TIP PRIMING

- The target surgical field must be properly prepared prior to either priming or applying BioGlue. BioGlue works best when the target surgical field is dry. A dry surgical field can be described as a field that does not retain with blood within 4-5 seconds after wiping dry with a surgical sponge.

CAUTION: Do not attempt to apply BioGlue to a field that is too wet. Application of BioGlue into a wet field may result in the failure of BioGlue to adhere.

- The residual air space must be removed prior to BioGlue application. To purge the air space, hold the Delivery Device with the syringe upright to ensure that the air bubbles in the solutions are located at the top of the syringe. Purging the air space can now be accomplished using two different methods:
 - Compress the Delivery Device trigger only until the solutions are even with the top of the syringe body. Once the residual air space has been removed, the syringe is ready for priming (refer to Step 3) and immediate use.
 - Compress the Delivery Device trigger until both solutions can be visibly seen in the base of the tip. The air space has now been removed, but this tip is now occluded with polymerized BioGlue and will need to be changed prior to priming (refer to Step 3) and application to the target site.

NOTE: Each syringe only needs to be purged of residual air space upon its initial use.

- Each applicator tip must be primed prior to BioGlue application. Priming ensures the BioGlue solutions are properly mixed. The surgeon should compress the Delivery Device trigger and expel a narrow ribbon of BioGlue approximately 3 cm long onto a sterile disposable surface (e.g., sponge, gauze, or towel).
- The surgeon should examine the material expelled during priming and ensure that it is of uniform light yellow to amber color and that it is free from air bubbles. If this material looks colorless or contains bubbles, repeat the prime as outlined in Step 2 until the device delivers a uniform liquid with no bubbles.

CAUTION: Avoid direct contact with material expelled during priming.

- When the applicator tip has been properly primed, proceed immediately to application.

CAUTION: BioGlue polymerizes very quickly. The surgeon must apply BioGlue immediately after priming. Pausing between priming and application can cause polymerization of BioGlue within the applicator tip. Should this occur, replace the obstructed tip with a new tip and repeat the steps for applicator tip priming. Do not continue to apply pressure to the trigger once the tip has occluded.

PRODUCT INFORMATION DISCLOSURE

Handling and storage of this device by the user as well as factors related to the patient, the patient's diagnosis, treatment, surgical procedures, and other matters beyond manufacturer's control may directly or indirectly affect this device and the results obtained from its use. This device should not be used except on the order of a physician.

DISCLAIMER OF IMPLIED WARRANTIES

CryoLife, Inc. disclaims all implied warranties with respect to BioGlue, including but not limited to, the implied warranty of merchantability and any warranty of fitness for a particular purpose. Prices, specifications, and availability are subject to changes without notice.

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**CryoLife[®], Inc.**
**CryoLife[®] Europa, Ltd.**