# BioGlue<sup>®</sup> | Surgical Adhesive

## Instructions for Use

BioGlue  
Kasutusjuhend  
Süstal

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- Bavaria J. E. *et al.* Advances in the treatment of acute type A dissection: an integrated approach. *Ann. Thorac. Surg.*, 2002, 74, S1848–52.
- Bavaria J. E. Acute aortic dissection: current perspectives. *J. Penn. Assoc. Thorac. Surg.*, 1999, 7–8.
- Coselli J. S. *et al.* Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J. Am. Coll. Surg.*, 2003, 197 (2), 243–253.
- Coselli J. S. *et al.* Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg.*, 2000, 5 (4), 259–76.
- Glock *et al.* Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue – CryoLife). *Laparoscopic Aorticliac Surgery for Occlusive Disease and Aneurysms*. 28. jaan. 2000.
- Hasse J. *et al.* Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J. Cardio-Thorac. Surg.*, 2001, 19 (1), 4–9.
- Kumar A. *et al.* Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (trükiis) *J. Clinical Neuroscience*.
- Kumar A. *et al.* Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J. Clinical Neuroscience*, 2003, 10 (1), 92–95.
- Lang H. *et al.* Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
- LeMaire S. A. *et al.* A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire S. A. *et al.* A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire S. A. *et al.* BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann. Thorac. Surg.*, 2002, 73 (5), 1500–5.
- LeMaire S. A. *et al.* Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire S. A. *et al.* Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J. *et al.* BioGlue surgical adhesive—an appraisal of its indications in cardiac surgery. *Ann. Thorac. Surg.*, 2002, 74, 432–7.
- Potaris K. *et al.* Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum*, 2003, 6 (5), 429–433.
- Potaris K. *et al.* Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor*, 2003, 9 (7), 179–83.
- Raanani E. *et al.* Use of “BioGlue” in aortic surgical repair. *Ann. Thorac. Surg.*, 2001, 72, 638–40.
- Sen A. *et al.* Cerebrospinal fluid leak rate after the use of BioGlue in translabirinthine vestibular schwannoma surgery: A prospective study. *Otol. Neurotol.*, 2006, 27, 102–5.
- Tansley P. *et al.* A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J. Thorac. Cardiovasc. Surg.*, 2006, 132, 105–12.
- Westaby S. *et al.* Acute type A dissection: conservative methods provide consistently low mortality. *Ann. Thorac. Surg.*, 2002, 73, 707–13.

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- Bavaria JE, *et al.* Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002;74:S1848-52.
- Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999;7-8.
- Coselli JS, *et al.* Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003;197(2):243-253.
- Coselli JS, *et al.* Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259-76.
- Glock *et al.* Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue – CryoLife). *Laparoscopic Aortolial Surgery for Occlusive Disease and Aneurysms*. Jan. 28, 2000.
- Hasse J, *et al.* Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001;19(1):4-9.
- Kumar A, *et al.* Evaluation of the use of BioGlue in neurosurgical procedures. 2003 (In Press) *J Clinical Neuroscience*.
- Kumar A, *et al.* Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92-95.
- Lang H, *et al.* Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). A pilot study. Presented at the XVIIIth European Association of Urology on 13 March 2003 (Abstract).
- LeMaire SA, *et al.* A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. Presented at the 34th Association of Academic Surgery 2000.
- LeMaire SA, *et al.* A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. Presented at Symposium for Cardiac and General Thoracic Surgery Teams. (Abstract) 2000.
- LeMaire SA, *et al.* BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002;73(5):1500-5.
- LeMaire SA, *et al.* Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- LeMaire SA, *et al.* Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. Presented at the 35th Association of Academic Surgery Meeting. (Abstract) 2001.
- Passage J, *et al.* BioGlue surgical adhesive—an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002;74:432-7.
- Potaris K, *et al.* Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Heart Surgery Forum* 2003;6(5):429-433.
- Potaris K, *et al.* Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. *Medical Science Monitor* 2003;9(7):179-83.
- Raanani E, *et al.* Use of "BioGlue" in aortic surgical repair. *Ann Thorac Surg* 2001;72:638-40.
- Sen A, *et al.* Cerebrospinal fluid leak rate after the use of BioGlue in translabyrinthine vestibular schwannoma surgery: A prospective study. *Otol Neurotol* 2006;27:102-5.
- Tansley P, *et al.* A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006;132:105-12.
- Westaby S, *et al.* Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002;73:707-13.

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- Bavaria JE, et al. Advances in the treatment of acute type A dissection: an integrated approach. *Ann Thorac Surg* 2002 m.;74:S1848–52.
- Bavaria JE. Acute aortic dissection: current perspectives. *J Penn. Assoc. Thorac Surg* 1999 m.;7–8.
- Coselli JS, et al. Prospective randomized study of a protein based tissue adhesive used as a hemostatic and structural adjunct in cardiac and vascular anastomotic repair procedures. *J Am Coll Surg* 2003 m.;197(2):243–253.
- Coselli JS, et al. Thoracic aortic anastomoses. *Operative Techniques in Thoracic and Cardiovascular Surg* 2000;5(4):259–76.
- Glock et al. Experimental techniques of aorto-prosthetic anastomosis by gluing (BioGlue – CryoLife). *Laparoscopic Aortoiliac Surgery for Occlusive Disease and Aneurysms*. 2000 m. sausio 28 d.
- Hasse J, et al. Experimental use of an albumin-glutaraldehyde tissue adhesive for sealing pulmonary parenchyma and bronchial anastomoses. *European J Cardio-Thorac Surg* 2001 m.;19(1):4–9.
- Kumar A, et al. Evaluation of the use of BioGlue in neurosurgical procedures. 2003 m. (spaudoje) *J Clinical Neuroscience*.
- Kumar A, et al. Reconstruction of the sellar floor using BioGlue following transphenoidal procedures. *J Clinical Neuroscience* 2003;10(1):92–95.
- Lang H, et al. Use of BioGlue in nephron sparing surgery for renal cell carcinoma (RCC). *Bandomasis tyrimas. Pristatyta XIII Europos urologų susivienijimo metu* 2003 m. kovo 13 d. (Santrauka).
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes acute phrenic nerve injury and diaphragmatic paralysis. *Pristatyta per Akademinės chirurgijos 34-ąjį susivienijimą* 2000 m.
- LeMaire SA, et al. A new surgical adhesive (BioGlue) causes immediate and prolonged bradycardia due to sinoatrial node degeneration. *Pristatyta širdies ir krūtinės chirurgijos simpoziumo dalyviams. (Santrauka)* 2000 m.
- LeMaire SA, et al. BioGlue surgical adhesive impairs aortic growth and causes anastomotic strictures. *Ann Thorac Surg* 2002 m.;73(5):1500–5.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects exposed nerves during the application of BioGlue Surgical Adhesive. *Pristatyta per Akademinės chirurgijos 35-ąjį susivienijimą. (Santrauka)* 2001 m.
- LeMaire SA, et al. Chlorhexidine Gluconate Gel protects the myocardium and sinoatrial node during the application of BioGlue Surgical Adhesive. *Pristatyta per Akademinės chirurgijos 35-ąjį susivienijimą. (Santrauka)* 2001 m.
- Passage J, et al. BioGlue surgical adhesive—an appraisal of its indications in cardiac surgery. *Ann Thorac Surg* 2002 m.;74:432–7.
- Potaris K, et al. Experience with an albumin-glutaraldehyde tissue adhesive in sealing air leaks after bullectomy. *Širdies chirurgijos forumas* 2003 m.;6(5):429–433.
- Potaris K, et al. Preliminary results with the use of an albumin-glutaraldehyde tissue adhesive in lung surgery. 2003 m. Medicinos mokslo monitorius;9(7):179–83.
- Raanani E, et al. Use of “BioGlue” in aortic surgical repair. *Ann Thorac Surg* 2001 m.;72:638–40.
- Sen A, et al. Cerebrospinal fluid leak rate after the use of BioGlue in translabrynthine vestibular schwannoma surgery: Būsimasis tyrimas. *Otol Neutrotol* 2006 m.;27:102–5.
- Tansley P, et al. A prospective, randomized, controlled trial of the effectiveness of BioGlue in treating alveolar air leaks. *J Thorac Cardiovasc Surg* 2006 m.;132:105–12.
- Westaby S, et al. Acute type A dissection: conservative methods provide consistently low mortality. *Ann Thorac Surg* 2002 m.;73:707–13.

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©2018 CryoLife, Inc. Kõik õigused kaitstud.
CryoLife, Inc.
CryoLife Europa, Ltd. • Bramley
Gamintojas: „CryoLife, Inc.“
©2018 „CryoLife, Inc“. Visos teisės sau-
„CryoLife, Inc.“ prekės ženklai. „Tyvek“ yra