CardioGenesis[®]

SoloGrip[®] III Handpiece

DELIVERY SYSTEM

Part of the CardioGenesis TMR Holmium:YAG Laser System

Catalog Number HP-SG3

INSTRUCTIONS FOR USE

- Caution:** Rx only. United States Federal law restricts this device to sale by or on the order of a physician (applies to U.S. sales only). United States federal law further restricts the use of this device to practitioners who have been trained in laser heart surgery including laser operation.
- Caution:** Use of this device is restricted to patients who have signed a procedure-specific consent form to ensure that the risks associated with TMR have been fully explained and understood (United States only).

IMPORTANT!

- The SoloGrip III Handpiece is a fragile surgical instrument. Use caution when removing contents from package and during use. Excessive stress or tension on the optical fiber contained in the handpiece may result in device damage or malfunction.
- Always place the laser console near the sterile field. Position the white fiber coupler near the operative site to minimize tension on the fiber when the handpiece is in use.
- Do not bend the optical fiber at sharp angles. If breaks or fractures appear in the optical fiber, immediately discontinue use and replace with a new handpiece.
- Each SoloGrip III Handpiece is supplied sterile for single use only. Do not reuse or resterilize. Any attempt to clean and resterilize can damage the fiberoptic.
- Inspect sealed, sterile package before opening. Product is sterile only in unopened, undamaged package. If package is opened or damaged, or if seal is broken, contents may not be sterile and may cause infection in the patient.

Sterilized using ethylene oxide

Do not use if package is damaged

Do not re-sterilize

Do not re-use

Caution

Consult instructions for use

Non-pyrogenic

Not made with natural rubber latex

Rx ONLY

1. DEVICE DESCRIPTION

The CardioGenesis TMR Holmium:YAG Laser System (CardioGenesis Laser System) is composed of a CardioGenesis Holmium:YAG laser console and a Handpiece Delivery System. The laser radiation emitted from this system has a wavelength of approximately 2.1 microns, which is in the mid-infrared (invisible) range of the electromagnetic spectrum. Water is the target absorber for this laser wavelength. This laser emits 200 microsecond laser radiation pulses at a 5 Hertz pulse repetition rate. The system limits nominal maximum average power to 7 watts (1.4 Joules/pulse). These pulses are not synchronized with the cardiac cycle, and there is no visible aiming beam.

The laser energy is delivered to the target tissue via the optical fiber with a diameter of approximately 1 mm. The optical fiber is integrated into the SoloGrip III Handpiece, which provides single-handed control of the optical fiber.

2. INDICATIONS FOR USE

Transmyocardial revascularization (TMR) with the CardioGenesis Laser System is indicated for treatment of stable patients with angina (Canadian Cardiovascular Society Class 4) refractory to medical treatment and secondary to objectively demonstrated coronary artery atherosclerosis and with a region of the myocardium with reversible ischemia not amenable to direct coronary revascularization.

3. CONTRAINDICATIONS

No contraindications known.

4. WARNING AND PRECAUTIONS

- Patients who were unstable (defined as unweanable from intravenous anti-anginal medications) experienced 11% (16/150) peri-operative mortality (surgery + 30 days) compared to 5% (7/132) in patients who did not require intravenous anti-anginal medications.
- In the Post Approval Study of 358 patients treated with either the SoloGrip III or PEARL 5.0 Handpiece, patients with ejection fractions $\leq 30\%$ were at higher risk of death compared to patients with ejection fractions $>30\%$.
- Do not treat the myocardium in the area of a left ventricular mural thrombus because of potential for the creation of emboli.

4.1 Explosions or Fire Hazard – Do not operate the laser in the presence of flammable substances, including gases, anesthetics, cleaning agents, combustible materials, or other volatile substances.

- Combustible or flammable materials (e.g., surgical drapes, gowns or gauze) in the surgical field may be ignited by laser radiation unless they are kept wet or moistened.
- Surround the surgical field with wet towels or wet gauze.
- Modify all other flammable materials to make them fire-retardant (e.g., flame resistant surgical drapes and gowns). Minimize oxygen exposure as oxygen increases the combustibility of materials exposed to the laser radiation.
- The SoloGrip III Handpiece may only be used with the CardioGenesis laser console.

4.2 Laser Radiation – The laser is classified as a Class IV laser product as defined in the Code of Federal Regulations, 21 CFR Section 1040.10(b)(11).

- Avoid exposure to laser radiation at all times during the installation and operation of this laser as direct or reflected radiation may damage skin or eyes.
- DO NOT LOOK DIRECTLY INTO THE Holmium:YAG BEAM as it can cause permanent ocular damage.
- Protect the patient's eyes by covering them with wet gauze or protective eyewear.
- All operating room personnel must wear protective eyewear with a minimum optical density of 3 at a wavelength of 2.1 μm when the laser is in use.

4.3 Physician Training

The CardioGenesis Laser System should only be used by properly trained surgeons.

Each SoloGrip III Handpiece is a fragile surgical instrument. Use caution when removing contents from packaging and during use. Excessive stress or tension on the optical fiber contained in the handpiece may result in device damage or malfunction.

4.4 Handling and Sterilization of the SoloGrip III Handpiece

- Inspect sealed, sterile package before opening. Product is sterile only in unopened, undamaged package. If package is opened or damaged, or if seal is broken, contents may not be sterile and may cause infection in the patient. Do not use product after "use by date" on the package.

- The SoloGrip III Handpiece is sterilized with EO gas and is intended **for single use only. Do not resterilize or reuse. Any attempt to clean and resterilize can damage the fiberoptic.** Resterilization may adversely affect proper mechanical function which could result in patient injury.
- Always place the laser console near the sterile field. Position the white fiber coupler near the operative site to minimize tension on the fiber when the handpiece is in use.
- Do not bend the optical fiber at sharp angles. If breaks or fractures appear in the optical fiber, discontinue use and replace with a new handpiece.

Biohazard – After use, handle and dispose of the SoloGrip III Handpiece as appropriate for a biohazard.

4.5 Precautions During TMR

- Avoid fluid loading the patient immediately prior to or during the TMR procedure, since fluid loading was associated with higher mortality in the clinical studies.
- The operator should pause 1/2-2 minutes between the creation of every 2-5 TMR channels since pauses may reduce the likelihood of ventricular arrhythmias following treatment.
- TMR treatment should be limited to the lower 2/3 of the left ventricle to avoid the mitral valve and disruption of the conduction system near the AV groove.
- A typical application of the device creates the lowest number of TMR channels yielding a density of 1 TMR channel/cm² in ischemic myocardium. The benefits of creating larger numbers of channels (> 40), either to treat broader areas of ischemia or to achieve a higher channel density, should be weighed against the potential risks of adverse events for each patient. In one IDE study of PEARL 5.0 & PEARL 8.0 Handpieces, 3 deaths occurred in 70 patients (4.3%). All 3 patients received 50 TMR channels. When all SoloGrip III Post Approval Study data were evaluated along with the PEARL 5.0 & PEARL 8.0 IDE data, patients with >40 TMR channels had a mortality of 6.0% compared to mortality of 1.7% for patients with ≤ 40 channels.
- The surgeon should have a defibrillator readily available throughout the surgical procedure.
- If the patient experiences ventricular fibrillation during the procedure, discontinue the procedure and treat the arrhythmia as appropriate.
- Patients receiving less reversible anti-coagulant therapy (e.g., Plavix®) should be changed to heparin at least 5 days prior to the procedure so that anti-coagulation can be more readily reversed.
- Patients should be weaned from antiplatelet therapy (e.g., Plavix) for at least 5 days prior to undergoing the TMR procedure.
- Anti-arrhythmic drugs should be considered for the immediate post-operative period to reduce the potential for arrhythmias.

5. ADVERSE EFFECTS

5.1 Observed Adverse Events

The randomized trial of TMR using the CardioGenesis Laser System with the SoloGrip III Handpiece versus medical management (MM) involved 275 patients who were followed for 204 patient-years.

There was one intraoperative death in the TMR group, which occurred in a patient who did not receive TMR – the patient developed ventricular fibrillation which could not be converted during the preparation of the heart for TMR. Within 30 days of TMR, five other patients died of cardiac causes and one died of pulmonary causes. In the MM group, two patients died within 30 days of enrollment in the study, both due to cardiac causes. During 12 month follow-up, an additional nine patients in the TMR group died (six due to cardiac causes, one each due to renal causes, multi-system organ failure and sudden death), and an additional five patients died in the MM group (all due to cardiac causes).

Adverse events were reviewed by an independent, masked Data and Safety Monitoring Board (**Table 1**).

5.2 Potential Adverse Events

Adverse events potentially associated with the use of TMR include (in alphabetical order):

- Accidental Laser Hit
- Acute Myocardial Infarction
- Arrhythmia
- Cerebrovascular Accident
- Conduction Pathway Injury
- Congestive Heart Failure
- Death
- Mitral Valve Damage
- Pulmonary Complications
- Unstable Angina

**Table 1: Percent and Number of Patients (n) with Adverse Events from the Pivotal Trial
TMR with SoloGrip III Handpiece versus Medical Management**

All patients in the randomized trial (n=275).

Recorded from 1996-1998.

Includes all adverse events, both related and unrelated to TMR, sorted alphabetically.

Adverse Event	TMR (n=132)				MM (n=143)			
	Early (0 to 30 days) % # events	Total (0 days to 1 year) % # events	Early (0 to 30 days) % # events	Total (0 days to 1 year) % # events				
Any Adverse Event	39% (51)	55% (72)	22% (31)	56% (80)				
Angina/Chest Pain Requiring Re-hospitalization	2.3% (3)	17% (22)	16% (23)	44% (63)				
Arrhythmia, Atrial	9.8% (13)	11% (14)	0.7% (1)	0.7% (1)				
Arrhythmia, Operative Ventricular Fibrillation (Op VF)	8.3% (11)	N/A	N/A	N/A				
Arrhythmia, Other Ventricular Arrhythmia	12% (16)	13% (17)	0% (0)	0% (0)				
Congestive Heart Failure	3.8% (5)	5.3% (7)	1.4% (2)	4.2% (6)				
Death (all causes)	5.3% (7)	13% ^a	1.6% (2)	8.6% ^a				
Dyspnea	0% (0)	0% (0)	1.4% (2)	8.4% (12)				
Hypotension	9.8% (13)	11% (14)	0% (0)	0% (0)				
Myocardial Infarction								
Q Wave MI	0.8% (1)	1.7% ^a	0.8% (1)	3.8% ^a				
Non Q Wave MI	4.5% (6)	12% ^a	0.8% (1)	6.7% ^a				
Pleural Effusion	0% (0)	2.3% (3)	0% (0)	0% (0)				
Respiratory Insufficiency	3.0% (4)	3.0% (4)	0% (0)	0% (0)				
Systemic Infection	1.5% (2)	1.5% (2)	0% (0)	0% (0)				
Transfusion Required								
Due to Blood Loss From TMR	0% (0)	N/A	N/A	N/A				
Due to Other Reasons	1.5% ^b (2)	1.5% (2)	0% (0)	0% (0)				
Unstable Requiring I.V. Anti-Anginals	1.5% (2)	17% (22)	19% (27)	48% (68)				

The following events were reported only once in the patients treated with TMR: allergic reaction, grand mal seizure, hemothorax, cardiomyopathy, pericarditis, peripheral edema, pneumothorax, pulmonary embolus.

The following events were reported only once in the patients treated with MM: cardiogenic shock, dehydration, pneumonia.

^a Kaplan-Meier estimate.

^b One due to GI bleed and one due to pre-existing anemia.

Note: Some patients experienced more than one adverse event.

**Table 2: Principal Safety and Effectiveness Results at 12 m from the Pivotal Trial
TMR with SoloGrip III Handpiece versus Medical Management**

All patients in the Randomized Trial (n=275).
Recorded from 1996-1998.

	TMR (n=132)	MM (n=143)	Difference (TMR-MM) [CI] ^a
Angina Improvement ^b	76% (58/76)	32% (16/50)	44%* [28%, 60%]
Thallium Scan Results ^c (n=61)			
Mean ± SD Δ Extent Ischemia %	-0.9 ± 9.4	-0.6 ± 10.8	-0.3 [-5.0, 5.6]
Mean ± SD Δ Extent Rest Defects %	1.6 ± 12.5	2.2 ± 11.8	-0.6 [-5.9, 7.1]
Freedom From All Cause Mortality			
30 Day Survival	95%	98%	3.7% [-1%, 8%]
Survival	87% ^d	91% ^d	4.9% [-2.5%, 12.3]
Event Free Survival ^e	55% ^d	31% ^d	24%* [12%, 35%]
Freedom From Treatment Failure ^f	74% ^d	48% ^d	26%* [16%, 38%]
Freedom From Hospitalization for Cardiac Causes	61% ^d	33% ^d	28%* [17%, 39%]
Medication Use			
Decrease in Calcium Channel Blockers (% Pts)	56%	24%	32%* [14%, 50%]
Decrease in Beta Blockers (% Pts)	39%	17%	22%* [6%, 39%]
Decrease in Nitrates (% Pts)	39%	24%	15% [-2%, 31%]
Quality of Life (DASI Score) ^g	21 ± 14	12 ± 11	9* [3.1, 14.9]
Exercise Treadmill Tests			
Total Exercise Time (min.)	7.9 ± 4.5	6.2 ± 5.6	1.7 [-0.6, 4.0]
Total Workload (METS)	5.0 ± 0.7	3.9 ± 0.8	1.1* [0.0, 2.1]

*= $p < 0.05$ P value calculated using Fisher's exact test, 2-sided for proportions, Student's t-test, 2-sided for continuous variables, or log rank test for Kaplan-Meier estimates.

^a CI = 95% confidence interval by normal approximation.

^b Angina Improvement: Improvement in angina symptoms from baseline to 12 months by (2) Canadian Cardiovascular Society Classes in patients who were available at 12 month follow-up.

^c Thallium Scans: A negative value indicates an improvement in a parameter. A positive value indicates a worsening.

^d Kaplan-Meier estimate.

^e Event Free Survival: Freedom from death, Q-wave MI, hospitalization for cardiac causes, CABG or percutaneous intervention.

^f Treatment Failure: Death, Q-wave MI, 2 cardiac hospitalizations within 3 months, 3 cardiac hospitalizations within 1 year, or unweanable from IV anti-anginal medications for at least 48 hours after at least two attempts at weaning.

^g DASI: Duke Activity Status Index for quality of life. A higher score indicates a better quality of life.

6. CLINICAL STUDIES

TMR VERSUS MEDICAL MANAGEMENT (MM)

Purpose: The purpose of this study was to compare TMR to MM. TMR was performed using the SoloGrip III Handpiece. **Primary outcome measures** were angina improvement, mortality, event-free survival, treatment failure, changes in perfusion as measured by thallium scans, and frequency of cardiac rehospitalizations.

Design and Patients: This multi-center, prospective, randomized, controlled trial was conducted at 18 U.S. centers. In the original study plan, 160 patients were enrolled between March 1996 and February 1997; 74 to TMR and 86 to MM. All 160 patients have reached one-year follow-up. Between February 1997 and July 30, 1998, an additional 115 patients (58 TMR, 57 MM) were enrolled in the study. Thus, a total of 275 patients were enrolled in the study through July 30, 1998; 132 were randomized to TMR and 143 were randomized to MM. Baseline characteristics and cardiac risk factors were similar between the two groups. Prior to 12 months, 46 patients in the MM group met a *priori* defined treatment failure criteria, became unstable, were withdrawn from this study, and were enrolled in a separate study of TMR use in unstable patients, leaving 97 patients in the MM group.

Methods: TMR was performed through a lateral thoracotomy. The handpiece, containing an optical fiber assembly, was applied to the surface of the left ventricle. As the fiber was advanced through the wall, pulsed laser energy was delivered to the myocardium until the optical fiber tip reached the left ventricle cavity. TMR channels were placed approximately 1 cm from each other. The number of channels created ranged from 16-87 (mean 39), using a mean energy of 1.4 Joules per pulse (range 1.2-1.6 J/pulse) and a mean of 14 pulses/channel.

Results: Table 2 lists the principal safety and effectiveness results. There were statistically significant differences in angina improvement and 12-month survival (event-free, freedom from treatment failure, and freedom from cardiac rehospitalization). There was no apparent difference in perfusion as measured by thallium scans (**Table 2**).

Kaplan-Meier estimates for mortality at 12 months were similar between the two groups: 87% for TMR treated patients and 91% for MM patients.

In the TMR group, five of 23 (22%) patients treated prior to July 1996 died within 30 days of the procedure. Investigators attributed these deaths to "fluid loading" patients in preparation for the TMR procedure, manipulation of the heart, and not pausing between creation of channels. These practices were altered in June 1996.

From July 1996 to completion of enrollment in July 1998, an additional 109 patients received TMR in the study. In this group, 30-day mortality was 1.8% (2/109).

TMR Post Approval Study

Purpose: The purpose of the study was to further define the disease characteristics of patients undergoing TMR and early (30-day) postoperative mortality predictors (risk factors). TMR was performed using the SoloGrip III Handpiece and the PEARL 5.0 Handpiece.

Design and Patients: This multi-center, single-arm, non-randomized post approval study was conducted at 18 US centers between March 1999 and February 2010. A total of 358 patients undergoing TMR therapy [286 retrospective (after PMA approval but before post approval study approval) and 72 prospective patients] were enrolled in the study. All patients at the time of the procedure were in CCS angina class IV and had an ejection fraction $\geq 25\%$ (mean=48.1 \pm 11.1). Endpoints for the study were early all-cause mortality and early major adverse cardiovascular events (MACE) including (cardiac related deaths, cerebrovascular accidents, Q-wave and non-Q wave myocardial infarctions (MI), congestive heart failure and serious arrhythmias). These data were compared to the pre-approval TMR data (n=132). Operator experience was characterized by assessing the relationship between the number of TMR procedures and patient mortality and MACE.

Methods: Medical history and procedure information were collected on all patients. Poolability of retrospective and prospective patients was assured through the enrollment procedures and similarities between demographics and medical histories. Patients were followed in the hospital and after discharge until 30 calendar days following the TMR procedure. During this 30-day period any mortalities or MACE were collected.

Results: Thirty day follow-up was complete for all patients (100%). There were 8 deaths (2.2%) reported within 30 days of TMR; 6 were cardiac related and 2 were non-cardiac related. There was no significant difference ($p=0.66$) between early mortality for retrospective (2.1%) and prospective (2.8%) patients. However, early mortality was statistically reduced ($p=0.0033$) in the post-approval patients (2.2%) compared with the pre-approval patients (5.3%). Univariate analysis of the 358 patient cohort indicated patients with ejection fractions $\leq 30\%$ (n=27) had a statistically significant increase in mortality compared to those patients with ejection fractions $>30\%$ (n=331; 11.1% vs. 1.5%, $p=0.0167$).

There were 43 MACE that occurred in 33 patients, resulting in a 30-day event rate of 9.2%. Serious arrhythmia accounted for the highest proportion of MACE (n=25, 58%), followed by congestive heart failure (n=8, 19%), cardiac related death (n=6, 14%), non-Q wave MI (n=2, 5%), Q-wave MI (n=1, 2%) and cerebral vascular accident (n=1, 2%). There was not a statistical difference in the 30-day MACE rate between prospectively and retrospectively enrolled patients (p=0.82). The number of TMR channels was the only statistically significant predictor of MACE in the study. Patients with ≥ 40 channels (n=83) had a statistically higher rate of MACE compared to those with < 40 channels (n=273; 18.7% vs. 7.0%, p=0.0095).

The surgeon's first two TMR patients (n=49) and the site's first three TMR patients (n=57) were considered as the learning curve. There was a slight but not statistically significant increase in early mortality during the surgeon's learning curve (4.1% vs. 1.9%, p=0.3010). However, the number of patients with at least one MACE during the surgeon's learning curve was statistically higher than the later period (18.4% vs. 7.8%, p=0.0293). The most notable difference was in arrhythmias with 16.3% occurrence in the learning curve compared to 5.2% for later patients (p=0.0091). There was not a statistical difference in MACE for the site's learning curve compared to later patients.

Strengths of the Post Approval Study included 100% follow-up of the enrolled patients and a large volume of centers demonstrating patients. Study weaknesses include the slow enrollment rate and use of retrospective patients.

7. PATIENT SELECTION AND TREATMENT

7.1 Specific Patient Populations

The safety and effectiveness of the CardioGenesis Laser System has not been established for the following specific populations:

- patients under the age of 18;
- patients who are pregnant or undergoing labor and delivery;
- nursing mothers;
- patients suffering from active hepatic disease, renal failure, cancer or major infection;
- patients with a left ventricular ejection fraction less than 25%;
- patients with mechanical heart valves;
- patients with angina CCS class III or better;
- patients with myocardial ischemia limited to the right ventricular wall.
- patients who have experienced a Q- wave myocardial infarction within the past 3 weeks.
- patients who have experienced a non-Q- wave myocardial infarction within the past 2 weeks.

8. PATIENT COUNSELING INFORMATION

This device is restricted to use in patients who sign an informed consent to ensure that the risks associated with TMR have been fully explained to, and understood by, the patient (United States only).

Patients should be advised that any reduction of angina may occur gradually, that they should continue on their anti-anginal medications, and that the need for these medications will be re-evaluated at subsequent visits.

Patients should be advised of the risks of the procedure including the possibility of:

- recurrence of angina;
- progression of myocardial ischemia;
- worsening heart failure;
- cardiac arrhythmia;
- death.

9. CONFORMANCE TO STANDARDS

The CardioGenesis Laser System has been tested to and conforms with the requirements of the following domestic and international standards:

- IEC 60601-1-1:2000 Medical Electrical Equipment, General Requirements for Safety;
- IEC 60601-1-2:2001 Medical Electric Equipment, General Requirements for Safety, Collateral; Electromagnetic Compatibility;
- IEC 60601-2-22:1995 Medical Electrical Equipment, Particular Requirements for the Safety of Diagnostic and Therapeutic Laser Equipment;
- 21 CFR 1040.10: U.S. Laser Product Performance Standards;
- EN 55011:1998 + A1:1999 (Radiated and Conducted Emissions);
- IEC 60825-1:1993 + A1:1997 + A2:2001 (Safety of Laser Products: Equipment Classification, Requirements, and User's Guide).

10. HOW SUPPLIED

10.1 Packaging

The CardioGenesis Laser System consists of the Holmium:YAG laser console and disposable Handpiece Delivery Systems.

The laser console is installed in the hospital by personnel authorized by CryoLife®.

The single use SoloGrip III Handpiece is supplied sterile. Do not use if package is damaged.

STERILE: Each SoloGrip III Handpiece is sterilized using EO and is supplied sterile for single use only. **Do not reuse or resterilize.** Resterilization may adversely affect proper mechanical function and could result in patient injury. Carefully inspect the sterile package for damage or defects prior to use.

PRECAUTION: Fragile surgical instrument. Use caution when removing contents from package. Excessive stress or tension on the optical fiber may result in device damage or malfunction.

10.2 Storage

The SoloGrip III Handpiece should be stored in a clean, dry environment under conditions that protect against extremes in temperature and humidity. To avoid device damage, do not stock other objects on device packaging. Proper stock rotation should be practiced.

10.3 Handling

The SoloGrip III Handpiece contains fragile fiberoptic elements. Always handle the instrument carefully. Place the laser console near the sterile field. Position the white fiber coupler near the operative site to minimize tension on the optical fiber when the handpiece is in use. Do not grip the optical fiber during testing or operation. Do not bend the optical fiber at sharp angles. If breaks or fractures appear in the optical fiber, immediately discontinue use and replace the device.

11. CLINICIAN USE INFORMATION

11.1 Patient Informed Consent

In addition to the standard surgical consent, a procedure-specific consent form is required to be signed by each patient prior to their procedure to ensure that the risks associated with TMR have been fully explained to the patient (United States only).

11.2 Device Operating Instructions

The SoloGrip III Handpiece has a plastic contoured handle with a thumbslide, which advances the optical fiber outward and retracts the optical fiber after formation of the TMR channel. Three depth stop settings allow the fiber to be advanced to three different lengths: 2.5 cm, 3.0 cm and 3.5 cm. Irrigation or a vacuum source can be used with the included stopcock valve and the SoloGrip III Handpiece as necessary.

The SoloGrip III Handpiece is designed for single-handed control of the optical fiber, which is incorporated into the device. The guideshaft of the SoloGrip III Handpiece is approximately 14 cm long and has an outside diameter of 3.8 mm. The distal guide shaft curves 90 degrees and the tip features a soft, flexible vacuum cup measuring 12 mm in diameter.

The distal portion of the SoloGrip III Handpiece handle features a rotation knob that allows for the manual rotation of the guide shaft. A vacuum source can be attached to the SoloGrip III Handpiece to enhance stabilization of the vacuum cup on the surface of the heart or irrigation can be channeled through the included stopcock valve to keep the epicardium moist.

11.3 Inspection Prior to Use

Carefully inspect the sterile package for damage or defects before opening. Do not use product after the Use By date of the package. If the integrity of the sterile package has been compromised or the package or product is defective, do not use the product and contact your CryoLife representative for return information.

PREOPERATIVE SET-UP

1. Test fire the laser prior to patient preparation.
 - a. Remove the CardioGenesis laser test fiber, supplied separately, from the package, remove the plastic end caps and inspect the fiber connectors before use.
 - b. Attach test fiber connector to laser aperture on the laser console. Fill a small non-metallic water basin with water approximately half full. Note: usage of too little water or a metal basin could result in damage to the test fiber or injury to the user.
 - c. Set laser to the calibration factor listed on disposable package label. Put on appropriate laser protective eyewear.

- d. Place the distal end of the test fiber so it just touches the surface of the water and activate the footswitch. Laser energy from the test fiber should cause a snapping sound.
 - e. Disconnect test fiber and dry the distal end of fiber. Replace both end caps and store the test fiber in its original packaging.
 - f. For more information on the test fiber, refer to the laser console Instructions For Use.
2. The circulating nurse should peel open the SoloGrip III Handpiece outer pouch to expose the inner sealed tray to the scrub nurse.
 3. The scrub nurse should peel the Tyvek® lid from the tray and pass off the fiberoptic connector to the circulating nurse for connection to the laser console.


Figure 1.

4. The circulating nurse should remove the plastic end cap and attach the fiberoptic connector to the laser console.
5. If vacuum is to be used, connect a sterile vacuum line to the stopcock valve at the distal end of the SoloGrip III Handpiece. Adjust vacuum to approximately 4-6 inches (102-152mm) Hg as desired.
6. If irrigation is to be used, connect a sterile normal saline solution to the stopcock valve at the distal end of the SoloGrip III Handpiece. Set the flow rate to approximately 1 drip per second.
7. The laser should be set up as explained in the laser console Instructions For Use.
8. Functionally test fire the SoloGrip III Handpiece in a container of sterile water or saline prior to use to ensure proper delivery of laser energy through the optical fiber. Ensure the optical fiber is in contact with the surface of the water/saline to protect the fiber while testing. Ensure the optical fiber does not come into contact with the walls of the container.


Figure 2.


CLINICAL USE

1. Key points to remember when using the SoloGrip III Handpiece:

- This is a delicate surgical instrument containing very thin fiberoptics.
- After lasing, retract the optical fiber completely before removing the vacuum cup from the epicardial surface.
- Do not bend the optical fiber at sharp angles when advancing or retracting through the myocardium.
- During the procedure, be alert for any sparking along the length of the fiberoptic. Flashing light in the coupler during use is normal.
- Failure to frequently inspect the condition of the optical fiber could lead to adverse events.

a. Place the vacuum cup securely against the epicardial surface.

b. Applying side-loaded pressure against the optical fiber while removing from the myocardium can damage the fiberoptic components.

Figure 3.

2. Vacuum may be used to stabilize the vacuum cup of the SoloGrip III Handpiece against the epicardium. The handpiece may also be used without the vacuum capability.
3. Use the depth stop lever to pre-select the optical fiber advancement length (either 2.5 cm, 3.0 cm or 3.5 cm). Push the depth stop lever to the right to select the 2.5 cm optical fiber advancement length. Push the depth stop lever to the left to select the 3.0 cm optical fiber advancement length. Leave the depth stop lever in the center position to select the 3.5 cm optical fiber advancement length.
4. To advance the optical fiber, push forward on the thumbslide while activating the footswitch to deliver the laser energy. Energy should be delivered to the treatment area until the optical fiber penetrates the ventricular wall. Penetration of the ventricular wall should require no more than 10-20 pulses from the laser console.
5. Following placement of each laser channel, retract the optical fiber completely inside of the handpiece before re-positioning the vacuum cup.
6. Initially and after every 2-5 channels are created, place the laser in "Standby" mode, fully extend the fiberoptics, clean and inspect the optical fiber. An indication of fiber breakage may be a discoloration. If discoloration is detected or if any damage to the optical fiber such as kinking, burning, melting, or breaking of the optical fiber is observed, discontinue use and replace with a new SoloGrip III Handpiece.


Figure 4.

Caution: During this inspection it is recommended that the laser console is placed on "Standby" mode.

7. If vacuum is being used, quickly break the vacuum seal by lifting the vacuum cup straight up as soon as the optical fiber has been retracted.
8. It is recommended that the SoloGrip III Handpiece is flushed or irrigated (maximum flow of 50 ml/minute) after approximately every 10 channels or when the thumbslide no longer moves freely.
9. The rotation knob may be rotated clockwise or counterclockwise to adjust the direction of the guideshaft and optical fiber.

11.4 Operator Training

Federal law restricts the use of this device to practitioners who have been trained in laser heart surgery including laser system operation. Operator training for use of the CardioGenesis Laser System must include training in the use of the laser console and the SoloGrip III Handpiece, as well as appropriate clinical training.

The CardioGenesis Laser System should only be used by personnel who have successfully completed training in the proper use of the CardioGenesis Laser System, as authorized by CryoLife. Training includes didactic and hands-on training in the following areas:

- Patient selection considerations;
- TMR technique;
- Patient management considerations.

Laser Training

The American National Standards Institute offers the following standard of practice for the use of lasers in medicine and surgery:

ANSI Z136.3, "American National Standard for Safe Use of Lasers in Health Care Facilities" 2005.

ANSI Z136.1, "American National Standard for Safe Use of Lasers" 2007.

For further information regarding training, contact 1-800-741-7062 (United States and Canada) and +1-770-419-3355 (International).

11.5 Mechanism of Action

The mechanism(s) whereby TMR relieves angina is not known. In addition to possible contributions of placebo effect, current theories include:

- Increased perfusion of myocardium via the channels created;
- Increased collateralization via angiogenesis;
- Symptom reduction resulting from disruption of pain fiber function;
- Possible microinfarcts to the myocardium.

12. PATIENT MANUAL


The brochure entitled "Information for Patients Considering Transmyocardial Revascularization" provides general information to the potential patient regarding the risks and benefits associated with TMR treatment.

13. PRODUCT INFORMATION DISCLOSURE

CryoLife has exercised reasonable care in the manufacture of this device. CRYOLIFE EXCLUDES ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS, SINCE HANDLING, STORAGE, REUSE AND RESTERILIZATION OF THIS DEVICE BY USER, AS WELL AS FACTORS RELATING TO THE PATIENT, DIAGNOSIS, TREATMENT, AND OTHER MATTERS BEYOND CRYOLIFE'S DIRECT CONTROL AFFECT THIS DEVICE AND THE RESULTS OBTAINED WITH ITS USE.

CRYOLIFE SHALL NOT BE LIABLE FOR ANY INCIDENTAL, OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS DEVICE. CRYOLIFE NEITHER ASSUMES, NOR AUTHORIZES ANY OTHER PERSON TO ASSUME FOR IT, ANY OTHER ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THIS DEVICE.

14. SYMBOLS

 = Sterilized using ethylene oxide	 = Do not re-use	 = Do not use if package is damaged
 = Caution	 = Batch code	 = Consult instructions for use
 = Catalogue number	 = Use-by date	 = Do not re-sterilize
 = Serial number	 = Not made with natural rubber latex	 = Non-pyrogenic
RxONLY	 = Quantity	 = Calibration factor

CryoLife[®]
Life Restoring Technologies[®]

www.CryoLife.com

Manufacturer

CryoLife, Inc.
1655 Roberts Boulevard, NW
Kennesaw, Georgia 30144
United States

1-800-741-7062, United States and Canada
1-770-419-3355, International
1-770-590-3753, Fax

CryoLife, Life Restoring Technologies, the snowflake design, CardioGenesis, and SoloGrip are trademarks of CryoLife, Inc. All other trademarks used are the property of their respective trademark holders.

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© 2019 CryoLife, Inc. All rights reserved.
© 2019 CryoLife, Inc. Tous droits réservés.