

CryoArtery[®] | Aortoiliac Artery

For Replacing Infected Aortic Grafts



CryoLife[®]
Life Restoring Technologies[®]



Outline

CryoArtery[®] | Aortoiliac Artery


- Indications for Use
- Underlying Need for CryoArtery Aortoiliac
- Target Patients
- Clinical Outcomes & Cost Savings
- CryoArtery Aortoiliac Artery vs. Alternative Procedures
- Cryopreserved Aortoiliac Allografts vs. Extra-Anatomic Bypass
- Implant & Post-Op Considerations
- Additional Benefits
- Arterial Allograft Configurations
- Extensions: Allograft Configurations
- Wet Lab Resources
- CryoLife Difference



Indications for Use

CryoArtery® | Aortoiliac Artery

- Indications for Use
 - Cryopreserved vascular allografts are indicated for use as vascular grafts.
- Contraindications
 - No contraindications for use of CryoArtery Aortoiliac Arterial Allografts are known.




Underlying Need

CryoArtery® | Aortoiliac Artery




Aortic Aneurysm Patient
102,256 ('16)¹




Endograft
84,097 ('16)¹

or



Open Surgery
18,159 ('16)¹


=



Infection
(0.2% to 5%)
204 to 5,112¹⁻³




Treatment Options




CryoArtery

or




NAIS

or



Ax-bifem

or




Rifampin-Soaked
Dacron

Target Patients: Abdominal Aortic Infections


CryoArtery® | Aortoiliac
Artery

Infected Synthetic Grafts




The graft infected is typically
(woven polyester “Dacron”
or ePTFE)

Aorto-Enteric Fistula



A fistula is created between the
bowel and synthetic graft. This
causes sepsis/infection in
the blood stream.


Mycotic Aneurysms



An infected aneurysm (picture
above: an abdominal aorta
with a mycotic aneurysm)

Clinical Outcomes & Costs Savings

CryoArtery® | Aortoiliac Artery



- 96% Freedom from Infection at 5 years²
- “Lower rates of graft rupture, aneurysm formation, recurrent infection and limb loss than other alternatives”²
- Potential cost and time savings in OR and ICU^{2,4}

**The Natural Choice
for Infected Fields**

CryoArtery Aortoiliac Artery vs. Alternative Procedures

CryoArtery® | Aortoiliac Artery

	CryoArtery Aortoiliac Artery 2,5,6,24	NAIS* Procedure ^{7,10}	Axillofemoral Reconstruction 2,8,9,11-17	Rifampin-Soaked Synthetic Graft ^{4,9}
Graft Re-infection	0% - 9%	10% - 14%	0% - 25%	7% - 47%
Mortality (30 day)	0% - 17%	0% -20%	11% - 28%	8% - 18%
Mean OR Time	4 – 7 hours	5 – 12 hours*	6 – 10 hours*	NR
Mean Length of Hospital Stay	16.7 – 24 days	21 – 28 days	18 – 33 days	30 days

*Some of these cases may have been performed as staged procedures.



“[...]Cryopreserved aortoiliac allografts should be considered a first line treatment against primary aortic graft infections.” (n=220)²

2. Harlander-Locke, et al. JVS 2014.
4. Vogt, et al. J Thoracic Cardiovasc 1998.
5. Vardanian, et al. AM Surg 2009.
6. Zhou W, et al. Tex Heart Inst J 2006.
7. Ali, et al. J Vasc Surg 2009.
8. Liedenbaum, et al. Worl J Surg 2009.

9. Bandyk, et al J Surg Res 2001.
10. Claggett, et al. J Vasc Surg 1997.
11. O'Hara, et al. J Vasc Surg 1986.
12. Reilly, et al. J Vasc Surg 1987.
13. Yeager, et al. J Vasc Surg 1999.
14. Seeger, et al. J Vasc Surg 2000.

15. Hart, et al. Ann Vasc Surg 2005.
16. Schmitt, et al. J Vasc Surg (Torino). 1990.
17. Bandyk, et al. J Vasc Surg 2001.
24. Brown KE, et al J Vasc Surg 2009.

Cryopreserved Aortoiliac Allografts vs. Extra-Anatomic Bypass⁴


	Arterial Allograft 	Extra-anatomic Bypass 	p value
Surgery/ Infection Related Mortality	12%	32%	0.008
Reoperation	9%	45%	0.001
Infection Completely Eliminated	91%	53%	0.001
Time in ICU	Median: 1 day Range: 1- 42 days	Median: 11 days Range: 2-120 days	0.002
Duration of Hospital Stay	Median: 14 days Range: 7-150 days	Median: 30 days Range: 15-240 days	0.02
Cost	Median: \$58,000 Range: \$55,000-\$160,000	Median: \$392,000 Range: \$89,000-\$580,000	0.005

“[...]Costs were 40% lower in the group treated by allografts.”⁴

Implant Technique Considerations

CryoArtery® | Aortoiliac Artery

Below is information to assist in matching the patient's pressurized aorta to a healthy unpressurized/undistended CryoArtery Aortoiliac Artery (AI)



*From measured unpressurized/undistended diameter to 120mmHg, (39% distension)²⁰


Implant Technique Considerations (continued)

CryoArtery® | Aortoiliac Artery

Below are illustrations of implant techniques used by vascular surgeons to suture an unpressurized/undistended CryoArtery Aortoiliac Artery to the native aorta to achieve an optimal diameter.²⁹




- Original CryoArtery Aortoiliac Graft



- More space between suture bites on the native aorta



- Use a beveled anastomosis to accommodate the difference in diameters



- “Fish Mouth” the CryoArtery Aortoiliac Graft to match the diameter of the native aorta



- Cut across the renal arteries and use the extra tissue as a taper to the native aorta



- Taper the native aorta to reduce diameter

Implant Technique Considerations (continued)

CryoArtery® | Aortoiliac
Artery



- Check for leaks: distend with D5LR to check & repair leaks (before implantation)¹⁹
 - Through and through ligature of the side branches
 - >2cm CryoLife does NOT suture ligate




- Atraumatic vascular clamp or/ & bulldog is/are only what should be used to occlude allograft




- Appropriate Length of Conduit¹⁸
 - Avoid kinking the allograft by anticipating the total length the allograft once it is pressurized
 - Absolute tension-free anastomosis

Implant Technique Considerations (continued)

CryoArtery® | Aortoiliac
Artery



- Anastomotic heel enlargement
- Allograft aortic portion distends up to 39% at 120mmHg²⁰
 - Begin with aortic anastomosis first**




- Anastomotic reinforcement with allograft strips¹⁸

Implant Technique Considerations (continued)

CryoArtery® | Aortoiliac
Artery

- Implant the posterior side aortoiliac graft to the anterior position to monitor for leaks from the lumbar arteries²⁸




- End-to-end 90° anastomosis when suturing the aortic portion of the allograft to the aorta using a single non-absorbable running polypropylene suture²³
- Spatulated 45° anastomosis when suturing the allograft iliac branches to the autologous iliac arteries²³
- Omentum wrapping around anastomosis to prevent fistula formation²²

Post-Op Protocol Considerations

CryoArtery® | Aortoiliac Artery

- Antibiotics:
 - When bacteria was identified, patients received pathogen-specific postoperative antibiotics for a minimum of 6 weeks²⁴
- Follow-Up:
 - Annual evaluation: to include imaging study (CT, MRA, Angiography, Doppler US)²⁴




- Aggressive wound drainage¹⁸

Additional Benefits

CryoArtery[®] | Aortoiliac Artery





- Natural Suturability²⁵
- Compliance at the Anastomosis



- Natural Pulsatile Flow²⁶

Arterial Allograft Configurations

CryoArtery® | Aortoiliac Artery

Allograft Type	Diameter	Length	Catalog #
Aortoiliac Artery 	Aorta: 8 mm – 15 mm*	6 cm – 11+ cm	R010
	Iliac: 4 mm – 5+ mm	4 cm – 11+ cm	
Descending Thoracic Artery 	8 mm – 15mm^	6 cm – 11+ cm	A020

*Diameter distends up to 39% at 120mmHg²⁰

^Diameter distends up to 60% at 120mmHg²⁷

Extensions: Allograft Configurations

CryoVein® | Vascular Allograft

CryoArtery® | Vascular Allograft

CryoVein®

Allograft Type	Diameter*	Length	Catalog #
Femoral Vein	6 mm – 15 mm	10 cm – 30+ cm	V060



CryoArtery®

Allograft Type	Diameter*	Length	Catalog #
Femoral Artery	4 mm – 5+ mm	10 cm – 30+ cm	R020





*DIAMETER: label changes as of June 15, 2015

- Beginning Monday, June 15, 2015, the following CryoLife cryopreserved femoral allografts mentioned above have been packaged and labeled with only the outer diameter on the package label and Certificate of Assurance (Note: the internal diameter is not referenced).
- Prior to Monday, June 15, 2015, the allografts listed above were labeled with only the internal diameter. The internal diameter was determined by measuring the outer diameter of the allograft and subtracting the estimated thickness for the vessel wall.
- In the future, customers may receive allografts that are labeled with the internal diameter or outer diameter (depending on the date the allografts were labeled).
- The Certificate of Assurance is the suggested reference document to review prior to ordering any of the allografts mentioned above to verify that the tissues meet your expectations.

Wet Lab Resources

CryoArtery® | Aortoiliac
Artery

Aortoiliac Implant Simulator Used in CryoLife Wet Labs



The CryoLife Difference

Experience

- Over 30 years of expertise in cryopreserving allografts (founded: 1984)
- Over 84,000 CryoLife vascular allografts shipped

Data

- 49 published clinical papers

Service

- Direct Representatives
- Hands-on Wet Labs

Quality

- Certified to be compliant with the ISO Quality System for Tissue Processing & Distribution
- Polypropylene monofilament suture for ligations (which does not harbor infection)
- Packaging: may be submerged in liquid nitrogen*
- CryoFreezers: available to hospitals for allograft storage
- AATB Accredited



*Prior to exposing the allograft to liquid nitrogen, allow the allograft to equilibrate with the liquid nitrogen vapor for a minimum of two (2) hours. After the equilibration period, the allograft may be submerged in liquid nitrogen for long-term storage. See Unpacking and Repacking Instructions for full details.



Learn more at www.CryoLife.com

Surgical technique is at the discretion of the surgeon. Variations in technique will inevitably and appropriately occur when clinicians take into account the needs of the individual patients.



CryoLife[®]
Life Restoring Technologies[®]