# BioFoam®

Surgical Matrix

## Instructions for Use

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CE0088

RxONLY

# **BioFoam® Surgical Matrix** **Instructions for Use**

Read the *Instructions for Use* prior to using this product.

## **DEVICE DESCRIPTION**

BioFoam® Surgical Matrix (BioFoam) is the reaction product of a protein solution with a cross-linker solution, both containing foaming agents. BioFoam is provided in a double-barreled container that is pre-filled with these two solutions. Upon mixing of the two solutions in the applicator tip attached to the container, both the foaming and cross-linking reactions begin. The foamy liquid expelled from the tip is adhesive, and sets up within seconds as soft foam.

The protein solution consists of bovine serum albumin (BSA) and sodium bicarbonate in water. The cross-linker solution is glutaraldehyde and acetic acid in water. These solutions are filled into the barrels of the container at a 4:1 volume ratio, respectively, packaged into an inner foil pouch and an outer poly/Tyvek® pouch, and terminally sterilized by gamma-irradiation. The material is expelled by advancing the pistons within the barrels of the container utilizing a manual plunger.

The product characteristics associated with clinical effectiveness include the adherence of the product to the surface of the tissue (cross-linking) and expansion of the product which generates a mixed-cell foam. The foam provides a mechanical barrier to decrease blood flow and pores for the blood to enter, leading to cellular aggregation and enhanced hemostasis.

## **INDICATIONS FOR USE**

BioFoam is intended for use as an adjunct to hemostasis in cardiovascular surgery and on abdominal parenchymal tissues (liver and spleen) when control of bleeding by ligature or conventional methods is ineffective or impractical.

## **CONTRAINDICATIONS**

- Not for patients with a known sensitivity to materials of bovine origin.
- Not for patients with a known sensitivity to glutaraldehyde.
- Not for intravascular use.

## **WARNINGS**

- Do not use BioFoam as a substitute for standard methods of surgical repair, such as sutures or staples.
- Do not expose valve leaflets or intracardiac structures to BioFoam.
- Do not use to obliterate a false lumen in the repair of aortic dissection; BioFoam application into a false lumen may result in re-dissection or distortion of the true lumen necessitating reoperation.
- Do not allow BioFoam in either the uncured or polymerized form to contact circulating blood. BioFoam entering the circulation can result in local or embolic vascular obstruction.
- Avoid contact with skin or other tissue not intended for application.
- Minimize use of BioFoam in patients with abnormal calcium metabolism (e.g., chronic renal failure, hyperparathyroidism). Glutaraldehyde-treated tissue has an enhanced propensity for mineralization.
- Do not use BioFoam if staff is not adequately protected (e.g., wearing gloves, mask, protective clothing, and safety glasses). Unreacted glutaraldehyde may cause irritation to eye, nose, throat, or skin, induce respiratory distress, and cause local tissue necrosis. Prolonged exposure to unreacted glutaraldehyde may cause a central nervous system or cardiac pathology. If contact occurs, flush affected areas immediately

with water and seek medical attention.

- Do not use BioFoam in the presence of infection and use with caution in contaminated areas of the body.
- Avoid repeat exposure of BioFoam in the same patient in subsequent procedure(s). Hypersensitivity reactions are possible upon exposure to BioFoam. Sensitization has been observed in animals.
- BioFoam contains a material of animal origin that may be capable of transmitting infectious agents.

## **PRECAUTIONS**

- It is recommended that surgical gloves, sterile gauze pads/towels, and surgical instruments be maintained moist to minimize the potential for BioFoam inadvertently adhering to these surfaces.
- BioFoam syringe and applicator tips are for single patient use only. Do not re-sterilize.
- Use only BioFoam applicator tips with the BioFoam syringe. The use of alternative applicator tips may negatively impact product performance (i.e., product expansion and adherence) or prohibit the appropriate application of the product to the target area.
- Do not use if packages have been opened or damaged.
- Take care not to spill contents of the syringe.
- Do not compress the syringe plunger while attaching it to the syringe.
- Do not apply BioFoam in a surgical field that is wet. This may result in poor adherence.
- Avoid tissue contact with material expelled from applicator during priming.
- BioFoam polymerizes rapidly. Priming must occur quickly, followed immediately by the application of BioFoam. Pausing between priming and application can cause polymerization within the applicator tip.
- Polymerized BioFoam has space occupying properties. Caution should be used in the application of BioFoam to avoid compression of the application site and adjacent structures.
- Allow BioFoam to polymerize for 2 minutes without manipulation. Manipulation of BioFoam prior to polymerization may result in a failure to provide hemostasis.
- To prevent the entrance of BioFoam into the cardiovascular system, avoid any negative pressure during application and polymerization of BioFoam. For example, (a) clamp and depressurize vessels prior to applying BioFoam to targeted anastomoses, or (b) turn off left ventricular vents prior to the application of BioFoam.
- BioFoam may constrict and deform veins. Appropriate steps should be taken to maintain lumen patency and shape, such as use of cannulae when applying circumferentially to veins.
- Do not use blood saving devices when suctioning excess BioFoam from the surgical field.
- The safety and effectiveness of BioFoam on patients in an anticoagulated state or patients with a coagulopathy have not been studied or established.

## **ADVERSE INCIDENTS**

Knowledge of all possible complications of cardiovascular and abdominal parenchymal (liver and spleen) surgery is necessary for physicians performing these procedures. Complications specific to these types of surgeries may occur at any time during or after the procedure.

Complications specific to the adjunctive use of BioFoam during cardiovascular and/or abdominal parenchymal (liver and spleen) surgery may include, but are not limited to the following: failure of product to adhere to tissue, failure of product to expand, failure of the device to dispense, failure of

product to achieve hemostasis, application of adhesive to tissue not targeted for procedure, vessel obstruction, thrombosis and thromboembolism, inflammatory and immune response, allergic reaction, mineralization of tissue, local tissue necrosis, injury to normal vessels or tissue, and possible transmission of infectious agents from material of animal origin.

## **PACKAGING AND STORAGE**

The BioFoam syringe and applicator tips are supplied sterile for single-patient use only. Discard any unused material from opened or damaged product.

The BioFoam solutions are contained within a capped, double-chambered, sterile syringe. Polymerized BioFoam is non-pyrogenic. Store at or below 25°C, but do not freeze.

## **DIRECTIONS FOR USE**


Apply BioFoam prophylactically or after a leak is detected.

## **Device Preparation**

The BioFoam delivery system consists of: solutions syringe, syringe plunger, and applicator tips.

Inside the BioFoam Solutions Syringe box there are two separate pouches. One contains the syringe and the syringe plunger, and one contains four applicator tips.


1. Remove the syringe, syringe plunger, and applicator tips from their packaging. While holding the syringe upright, tap the syringe until the air bubbles in the solutions rise to the top of the syringe.

**Figure 1**


**NOTE:** Continue to hold the syringe upright during the entire assembly of the delivery system to keep the bubbles toward the top of the syringe.

2. Remove an applicator tip from its packaging and inspect the collar portion of the tip to ensure that the pointer portion is directly adjacent to the larger port. If not, rotate the locking collar on the shaft until the pointer is adjacent to the larger port.

**Figure 2**

3. While firmly grasping the syringe, nose upward, turn the cartridge cap 90° counterclockwise and remove the cap by rocking it from side-to-side. Align the tip with the cartridge using the corresponding notches on each and place the tip on the syringe.

**Figure 3**

**CAUTION:** Take care not to spill solution from the syringe during assembly.

4. Lock the applicator tip in place by pushing the tip firmly toward the syringe and rotating the tip collar 90° clockwise.

Figure 4

5. While keeping the syringe upright, align the small and large barrels of the syringe with the corresponding syringe plunger heads and slide the plunger into the back of the syringe until resistance is felt. The syringe delivery device is now assembled.


Figure 5

**CAUTION:** Do not lay the assembled device on its side until all air has been purged (see next paragraph).

**CAUTION:** Before using BioFoam in the procedure, the syringe must be purged of the residual air space and the applicator tip must be primed. Refer to [Site Preparation, Syringe Air Space Removal and Applicator Tip Priming](#).

6. To remove occluded applicator tips, grasp the applicator tip collar, rotate the tip collar counterclockwise, and lift the tip off the syringe by rocking it side-to-side.

#### **Site Preparation, Syringe Air Space Removal and Applicator Tip Priming**

1. The target surgical field must be properly prepared prior to removal of the residual air space, priming, or applying BioFoam. BioFoam works best when the target surgical field is dry. A dry surgical field can be described as a field that does not retain with blood within 4-5 seconds after wiping dry with a surgical sponge. The surgeon should use methods such as digital pressure, clamping, etc. to obtain a temporary dry field prior to applying BioFoam.

**CAUTION:** Do not attempt to apply BioFoam to a field that is not dry. Application of BioFoam into a field that is not dry may result in the failure of BioFoam to adhere, resulting in potential formation of hematoma beneath the device surface. Insufficient adhesion of BioFoam to the tissue surface is identified by observance of blood leaking from the BioFoam-tissue interface. If this occurs, remove the non-adhered BioFoam and reapply to a dry surgical field (i.e. field does not retain with blood and/or bile within 4-5 seconds after wiping dry with a surgical sponge).

**NOTE:** The presence of water, irrigation fluids, blood, serum, thrombus, or bile will interfere with the adherence of BioFoam to the tissue surface.

2. The residual syringe air space must be removed prior to BioFoam application. Again, it is important to hold the assembled syringe upright to ensure that the air bubbles in the solutions are located at the top of the syringe. Purging of the air space can now be accomplished using two different methods:

- a. Compress the plunger only until the solutions are even with the top of the

syringe body. Once the residual air space has been removed the syringe is ready for priming (refer to Step 3) and immediate use.

- b. Compress the plunger until both solutions can be visibly seen in the base of tip. The air space has now been removed, but this tip is now occluded with polymerized BioFoam and will need to be changed prior to priming (refer to Step 3) and application to the target site.

Figure 6

**NOTE:** Each syringe only needs to be purged of residual air space upon its initial use.

3. Each applicator tip must be primed prior to BioFoam application. Priming ensures the BioFoam solutions are properly mixed. The surgeon should compress the plunger and expel a narrow ribbon of BioFoam approximately 3 cm long onto a sterile disposable surface (e.g., sponge, gauze, or towel).
4. The surgeon should examine the material expelled during priming and ensure that it is of uniform color and consistency. If this material does not foam or looks colorless, repeat the prime as outlined in Step 2 until the device delivers a uniform product.

**CAUTION:** Avoid direct contact with material expelled during priming.

5. When the applicator tip has been properly primed, proceed immediately to application.

**CAUTION:** BioFoam polymerizes very quickly. The surgeon must apply BioFoam immediately after priming. Pausing between priming and application can cause polymerization of BioFoam within the applicator tip. Should this occur, replace the obstructed tip with a new tip and repeat the steps for applicator tip priming (refer to Steps 3 and 4). Do not continue to apply pressure to the plunger once the tip has occluded.

#### **General Techniques for the Use of BioFoam in Surgery**

**NOTE:** Before using BioFoam, surgeons should become familiar, by appropriate training, with the surgical techniques and variations of their specific procedures. The use of BioFoam should be practiced with the product prior to initial use in the surgical suite.

1. BioFoam works best when the target surgical field is dry. A dry surgical field can be described as a field that does not retain with blood within 4-5 seconds after wiping dry with a surgical sponge.

**NOTE:** In an *in vitro* study measuring the tensile adhesive strength of BioFoam, no statistical difference was observed when BioFoam was applied onto a collagen membrane equilibrated at 10°C and 20°C.

2. Tissues surrounding the target surgical site should be protected from the unintentional application of BioFoam. The most effective method of protection is to cover any non-target tissues with moist sterile gauze pads.
3. Apply an even coating of BioFoam to the target area in a slow and steady manner in a quantity sufficient to provide hemostasis.

- a. Application of BioFoam to the tissue surface should be limited to one layer thick, which is approximately 0.5 cm tall.

- b. During application, the applicator tip should be maintained at a distance of ≤ 6 mm from the tissue surface, but not touching the tissue.

**CAUTION:** Avoid contact of BioFoam with blood-saving devices, such as cell savers and pumps.

**NOTE:** Avoid contact of the applicator tip with unpolymerized BioFoam once delivered to prevent polymerization of BioFoam on the applicator tip.

4. If BioFoam is inadvertently applied to non-target tissues, allow the adhesive to polymerize completely. Then, using forceps, carefully remove the polymerized BioFoam from the unintended area.

**CAUTION:** Polymerized BioFoam has space occupying properties. Caution should be used in the application of BioFoam to avoid compression of application site or adjacent structures.

5. BioFoam does not require any clamping or compression in order to polymerize.

**CAUTION:** Allow BioFoam to polymerize for 2 minutes without manipulation. Manipulation of BioFoam prior to polymerization may result in a failure to provide hemostasis.

6. When BioFoam has completely polymerized, the surgeon may trim away any excess material or irregular edges with scissors and forceps.

#### **General Techniques for the Use of BioFoam in Abdominal Parenchymal Surgery**

1. BioFoam should be used as an adjunct to standard methods of surgical repair in abdominal parenchymal surgery. Vessels ≥ 2 mm should be ligated prior to the application of BioFoam.
2. BioFoam should be applied from edge to edge of the resected tissue surface. BioFoam will expand approximately 0.5 cm beyond the edge.

#### **General Techniques for the Use of BioFoam in Cardiovascular Surgery**

1. BioFoam should be used as an adjunct to standard methods of cardiovascular surgical repair in tissue approximations.
2. BioFoam should be applied across the entire length of a suture line. BioFoam will expand approximately 0.5 cm beyond the suture line.

**CAUTION:** BioFoam may constrict and deform veins. Appropriate steps should be taken to maintain lumen patency and shape, such as use of cannulae when applying circumferentially to veins.

**CAUTION:** Do not expose valve leaflets or intracardiac structures to BioFoam.

**CAUTION:** Do not use to obliterate a false lumen in the repair of aortic dissection; BioFoam application into a false lumen may result in re-dissection or distortion of the true lumen necessitating reoperation.

**CAUTION:** To prevent the entrance of BioFoam into the cardiovascular system, avoid any negative pressure during application and polymerization of BioFoam. For example, (a) clamp and depressurize vessels prior to applying BioFoam to targeted anastomoses, or (b) turn off left ventricular vents prior to the application of BioFoam.

**RxONLY** = Prescription Device

= Latex Free

= Non-pyrogenic

= Do not use if package is damaged

In the United States of America Federal law restricts a prescription device to sale by or on the order of a physician.

In all other regions of the world, this device should not be used except on the order of a physician.

### **PRODUCT INFORMATION DISCLOSURE**

Handling and storage of this device by the user, as well as factors related to the patient, the patient's diagnosis, treatment, surgical procedures, and other matters beyond manufacturer's control, may directly or indirectly affect this device and the results obtained from its use. This device should not be used except on the order of a physician.

### **DISCLAIMER OF WARRANTIES: LIMITS OF LIABILITY**

CRYOLIFE DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES WITH RESPECT TO THIS SURGICAL MATRIX, INCLUDING BUT NOT LIMITED TO THE EXPRESS AND IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL CRYOLIFE BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. In the event that such disclaimer is found invalid or unenforceable for any reason: (i) any action for breach of warranty must be commenced within one year after any such claim or cause of action accrued and (ii) the remedy for any such breach is limited to the replacement of the product.

Prices, specifications, and availability are subject to change without notice.

Representative for vigilance reports (except for European Union countries):

CryoLife, Inc.

1655 Roberts Boulevard, NW

Kennesaw, Georgia 30144 · USA

Phone: 1-800-438-8285 (United States and Canada)

Phone: +1-770-419-3355 (United States and International)

Fax: +1-770-590-3753

Website: [www.cryolife.com](http://www.cryolife.com)

E-mail: [usbiofoam@cryolife.com](mailto:usbiofoam@cryolife.com)

CryoLife Europa, Ltd.

Bramley House • The Guildway • Old

Portsmouth Road

Guildford • Surrey GU3 1LR • United Kingdom

Phone: +44 (0) 1483 441030

Fax: +44 (0) 1483 452860

E-mail: [Europa@cryolife.com](mailto:Europa@cryolife.com)

Manufactured by: CryoLife, Inc.

**CryoLife**®

Life Restoring Technologies®

© 2012 CryoLife, Inc. All rights reserved. The snowflake design, CryoLife, Life Restoring Technologies and BioFoam are trademarks of CryoLife, Inc. Tyvek is a registered trademark of E.I. du Pont de Nemours and Company.

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CryoLife Europa Ltd. • Bramley