

The On-X[®] Experience

Medical Carbon Research Institute, LLC

September 2000/Volume 2/Issue 1

From the Editor

Robert More

Vice President Research and Development

Medical Carbon Research Institute (MCRI)

This issue of "The On-X[®] Experience" provides an update of clinical implantations of the On-X Valve, a report on the regulatory status, a comparison of linearized complication rates, a listing of recent publications and presentations and ancillary studies on the On-X Valve and On-X Carbon.

On-X Valve Update

The first clinical implant of an On-X[®] Prosthetic Heart Valve occurred on September 12, 1996. Since that time, approximately 3500 On-X valves have been implanted worldwide. Total valve implant experience now exceeds 3300 patient years.

The number of implanting centers continues to grow and is now more than 200 centers. The On-X valve has now been implanted in 33 countries. Recently, centers in Belgium, Slovakia, Croatia, Iraq and Iran have begun implanting the On-X valve. Implants in the US clinical centers continue under the IDE granted by the US FDA in December 1998.

Regulatory Status

The CE mark was obtained for the On-X valve from the TÜV in Munich on July 24, 1998. Regulatory approval to allow implantation of the On-X valve has been obtained in other countries as required. Applications are pending in Japan, Canada and the United States.

MCRI's interaction with the US FDA continues to proceed favorably. On Sept. 16, 1999, the FDA accepted the first module of MCRI's modular Pre-Market Approval (PMA) Application. The final module

was submitted August 30, 2000. The FDA staff will require some time to complete the application review. Clinical study data submitted to the FDA are summarized as follows.

Two Year Follow-Up Data

Patient population and surgery

From September 1996 to May 2000, 301 patients had isolated valve replacement at 11 European centers under a standardized protocol. Isolated aortic valve replacement (AVR) occurred in 184 (61%) patients and isolated mitral valve replacement (MVR) in 117 (39%). There were 21 double valve implants to give a total population of 322 patients.

In AVR, 66% of patients were male, while in MVR 46% were male. The mean age at implant was 60.2 ± 8.4 years in AVR (range 20-80) and 60.0 ± 10.2 years in MVR (range 21-76).

Preoperatively, approximately half of the AVR patients were in NYHA Class II, while 62% of the MVR patients were in Class III indicating that the mitral population was more severely ill to start. Additionally, only 2.2% of AVR patients had previous cardiac surgery, while 17.1% of MVR patients had previous cardiac surgery.

The most common disease etiology was calcific degeneration in AVR (50.0%) and rheumatic heart disease in MVR (38.5%). Other common etiologies were degenerative: 27.7% AVR and 27.4% MVR and endocarditis: 4.4% AVR and 11.1% MVR. The valve lesion leading to replacement was 46.7% stenosis, 21.2% regurgitation and 32.1% mixed in AVR and 17.1% stenosis, 51.3% regurgitation and 31.6% mixed in MVR. Typical of these populations, 87.5% of AVR patients were in sinus rhythm at surgery, while 53.9% of MVR patients were in atrial fibrillation. In both AVR

and MVR, all demographic variables examined indicated a usual adult mechanical valve replacement population.

Follow-up and data analysis

Follow-up in the study was obtained entirely through office visits. In accordance with the protocol, patients were maintained on anticoagulant therapy with a target therapeutic international normalized ratio (INR) range of 2.5 to 3.5 for AVR. For MVR, the INR target range was 3.0 to 4.5. The target range of anticoagulant therapy was maintained unless a contraindication to such therapy arose. A review of the INR measurements showed a median INR of 2.49 AVR. Mitral INR measurements showed a median INR of 2.68. For both the aortic and mitral measurements one-third were below 2.0. Thus, the On-X valve is already performing well in a population with low-end anticoagulant levels without apparent increase in TE or thrombosis rates.

Total follow-up in the study was 405.6 patient years (pt-yrs) AVR and 219.5 pt-yrs MVR. The average follow-up for all patients was about 24 months, 48 patients have been followed 3 years.

Morbid events were identified and classified in accordance with the AATS/STS guidelines for cardiac valve operations. Linearized rates for morbid events were calculated as percent per patient-year of follow-up (or number per 100 years of follow-up). Morbid event rates were compared to rates from the literature for other mechanical valves over comparable follow-up periods.

Results

Morbid events related to thrombosis and bleeding for isolated replacements are summarized in Table 1. The On-X valve has lower rates of thrombosis, thromboembolism and hemorrhage compared to literature reports on other mechanical valves.

Medical outcome

Outcome of valve replacement was determined by examining the preoperative and postoperative NYHA classification of each patient. Cross-tabulation at the 2 year interval shows that 75.8% of AVR patients and 70.6% of MVR patients improved one or more classes, 21.2% of AVR patients and 17.6% of MVR patients stayed the same and 3% AVR and 3.9% MVR

showed a higher classification.

Quality of Life

Quality of life surveys demonstrated that both AVR and MVR patients with preoperative serious complicated disease states were restored to a minor, uncomplicated chronic disease state, especially with regard to physical function.

On-X Valve, Bioprosthetic Stented and Stentless Valve Hemodynamics

The purported superior hemodynamics of some stented and stentless bioprosthetic valves has been used to rationalize implantation of these valves rather than mechanical valves. However, a detailed comparison of hemodynamic results for bioprosthetic valves (from the literature) to those of the On-X valve strongly suggests that this prima facie assumption of superior bioprosthetic hemodynamics is not true. In fact, the hemodynamic performance of the On-X valve was superior, particularly in the small aortic root sizes. Comparison results are given in Tables 2 through 4. Hemodynamic study citations used in this comparison were all performed in resting patients using the full Bernoulli equation form and the velocity-time integral form of the continuity equation.

Recent On-X Publications/Presentations **Performance of the On-X[®] Valve in the Small Aortic Root: A Multicenter Trial**

A. Haverich, and the Multicenter Investigators
Medizinische Hochschule Hannover, Hannover, Germany

Presented to the Asian Society for Cardiovascular Surgery, Fukuoka, Japan, September 6 – 8, 2000

A clinical trial of the performance of the On-X[®] Prosthetic Heart Valve, particularly in the small aortic root, was conducted at 11 centers in Europe. From September 1996 to January 2000, 184 aortic valves were implanted, including 17 19mm and 35 21mm valves. All patients were followed with a mean of 11 months. Hemodynamic values were measured at discharge and 1 year, and blood studies were done at each visit. Clinical outcome was evaluated by NYHA classification. Mean pressure gradient for the valve at dis-

charge was 12.6 ± 4.7 and 9.0 ± 5.3 mmHg for 19 and 21mm valves, while their effective orifice areas were 1.5 ± 0.2 and 1.8 ± 0.4 cm². Blood damage was low with LDH levels of 214 ± 46 - 196 ± 33 for size 19 and 212 ± 33 - 215 ± 33 for size 21 at 3-6 months and 1 year (upper normal = 250). NYHA improvement occurred in 75% of patients. The early performance of the On-X valve in the small aortic root is exceptional.

A Comparison of the Initial European and North American Experience With the On-X® Valve

A. Laczkovics, and the Multicenter Investigators
Klinikum Bergmannsheil, Ruhr University-Bochum,
Bochum, Germany

Presented to The International Society of Cardiothoracic Surgeons, Vancouver, August 13-16, 2000

A multicenter clinical investigation of the On-X® Prosthetic Heart Valve was conducted in Europe (EU) and an identical trial was started in North America (NA). Through December 1999, in the EU trial 301 patients received isolated valve implants at 11 centers and 61 isolated valve patients were entered in the NA trial at 9 centers. Follow-up averages about 2 years in EU and 1 year in NA. Demographically, the patients are similar: mean age 60.2 EU and 57.1 NA; 56% male EU and 50% male NA; primary etiology calcific degeneration for aortic valves and rheumatic or degenerative for mitral valves in both studies; valve lesion and preoperative New York Heart Association classification were also not different between studies. Follow-up is limited in the NA study, so comparison of adverse event rates is not feasible, but the rates are low in each study, for example, in the EU thromboembolism 2.2%/ptyr, paravalvular leak 0.7, and endocarditis 1.4. The principle early comparisons between the studies come from postoperative echocardiographic hemodynamics and blood studies for hemolysis. Aortic effective orifice area in EU ranged from 1.5 to 2.7 cm² for 19 to 25 mm valves, and the results are similar in NA at 1.4 to 3.4 cm². Postoperative lactate dehydrogenase at 3-6 months was 225 aortic and 268 mitral EU, and 212 aortic and 272 mitral NA (upper normal 250 I.U.). These two studies show similar patient populations, have excellent agreement in results and demonstrate the satisfactory early performance of the valve.

1. Chambers J, et al. *Journal Heart Valve Disease* 7:569-573,1998.S.
2. Fraund, et al. *Thoracic Cardiovascular Surgeon* 46:293-297,1998
3. Birnbaum D, et al; *Journal Heart Valve Disease*; 9:142-145,2000

Examination of Hemolytic Potential with the On-X Prosthetic Heart Valve

Birnbaum D, Laczkovics A, Heidt M, Oelert H, Laufer G, Greve H, Pomar J, Mohr F, Haverich A, Regensburger D.

The Journal of Heart Valve Disease 2000;9:142-145.

Comparison of the Classical and Modified Forms of the Continuity Equation in the Aortic Position

Chambers J, Ely J, members of the On-X Prosthetic Heart Valve Trial.

The Journal of Heart Valve Disease 2000;9:299-302.

In Vitro Hydrodynamics of Four Bileaflet Valves In Mitral Position

Feng ZG, Umezu M, Fujimoto T, Tsukahara T, Nurishi M, Kawaguchi D.

Asian Cardiovas Thorac Ann 2000;8:3-10.

Comparison of On-X and SJM HP bileaflet aortic valves

Walther T, Falk V, Tigges R, Krüger M, Langebartels G, Diegeler A., Autschbach R, Mohr F.

The Journal of Heart Valve Disease 2000;9:403-407.

Unalloyed Pyrolytic Carbon for Implanted Mechanical Heart valves

Ma Ling, Sines G.

The Journal of Heart Valve Disease 1999;8:578-585.

Assessment of Perceived Mechanical Heart Valve Sound Level in Patients

Nygaard H, Johansen P, Riis C, Hasenkam JM, Paulsen PK.

The Journal of Heart Valve Disease 1999;8:655-661.

Table I
Linearized Rate (% / pt.yr.) Comparison

Morbid Event	On-X:		SJM:		CMI:		Med-Hall:	
	Aortic	Mitral	Aortic	Mitral	Aortic	Mitral	Aortic	Mitral
Thromboembolism	1.73	1.82	1.71	1.86	1.11	2.05	1.50	2.00
Thrombosis	0.00	0.00	0.19	0.14	0.03	0.47	0.20	0.40
Hemorrhage	0.74	0.00	2.00	1.34	1.82	1.92	0.80	1.90
TOTALS	2.47	1.82	3.90	3.34	2.96	4.44	2.50	4.30

Rate Comparison References

- Edmunds Jr LH, Clark RE, Cohn LH, et al. Guidelines for reporting morbidity and mortality after cardiac valvular operations. *Ann Thoracic Surg* 1996;62:932-35
- Grunkemeier GL, Starr A, Rahimtoola, SH. Prosthetic heart valve performance: long-term follow-up. *Current Problems in Cardiol* 1992;17(6):333-406
- Jamieson WRE, Edwards FH, Schwartz M, et al. Risk stratification for cardiac valve replacement. *National Cardiac Surgery Database. Ann Thorac Surg* 1999;67:943-51
- Aagaard J, Hansen CN, Tingleff J, Rygg I. Seven-and-a-half years clinical experience with the CarboMedics Prosthetic Heart Valve. *J Heart Valve Disease* 1995;4:628-33
- Aoyagi S, Oryoji A, Nishi Y, et al. Long-term results of valve replacement with the St. Jude Medical valve. *J Thorac Cardiovasc Surg* 1994;108:1021-29
- Baudet EM, Puel V, McBride JT, et al. Surgery for acquired heart disease, long-term results of valve replacement with the St. Jude Medical prosthesis. *J Thorac Cardiovasc Surg* 1995;109:858-70
- Bech-Hanssen O, Wallentin I, Larsson S, et al. Reference Doppler echocardiographic values for St. Jude Medical, Omnicarbon, and Biocor prosthetic valves in the aortic position. *J Am Soc Echocardiogr* 1998;11:466-77
- Bernal JM, Rabasa JM, Gutierrez-Garcia F, et al. The CarboMedics valve: experience with 1,049 implants. *Ann Thorac Surg* 1998;65:137-43
- Bitar JN, Lechin ME, Salazar G, Zoghbi WA. Doppler echocardiographic assessment with the continuity equation of St. Jude Medical mechanical prostheses in the mitral valve position. *Am J Cardiol* 1995;76:287-93
- Copeland, III JG. An international experience with the CarboMedics Prosthetic Heart Valve. *J Heart Valve Disease* 1995;4:56-62
- Copeland JG. The CarboMedics Prosthetic Heart Valve in the mitral position: results of the multicenter international trial. *J Card Surg* 1997;12:205-09
- Craver J. CarboMedics Prosthetic Heart Valve. *Eur J Cardio-thoracic Surg* 1999;15(Suppl D):S3-S11
- Debetaz LF, Ruchat P, Humi M, et al. St. Jude Medical Valve Prosthesis: an analysis of long-term outcome and prognostic factors. *J Thorac Cardiovasc Surg* 1997;113:134-48
- De Luca L, Vitale N, Giannolo B, et al. *J Thorac Cardiovasc Surg* 1993;106:1158-65
- Fiane AE, Geiran OR, Svennevig JL. Up to eight years' follow-up of 997 patients receiving the CarboMedics Prosthetic Heart Valve. *Ann Thorac Surg* 1998;66:443-48
- Horstkotte D, Schulte H, Bircks W, et al. Unexpected findings concerning thromboembolic complications and anticoagulation after complete 10-year follow-up of patients with St. Jude Medical prostheses. *J Heart Valve Disease* 1993;2:291-301
- Ibrahim M, O'Kane H, Cleland J, et al. The St. Jude Medical prosthesis - a thirteen-year experience. *J Thorac Cardiovasc Surg* 1994;108:221-30
- Isomura T, Hisatomi K, Hirano A, et al. The St. Jude medical prosthesis in the mitral position. *Eur J Cardio-Thorac Surg* 1994;8:11-14
- Jamieson WRE, Miyagishima RT, Grunkemeier GL, et al. Bileaflet mechanical prostheses performance in mitral position. *Eur J Cardio-Thorac Surg* 1999;15:786-94
- Jegaden O, Eker A, Delahaye F, et al. Thromboembolic risk and late survival after mitral valve replacement with the St. Jude Medical valve. *Ann Thorac Surg* 1994;58:1721-28
- Khan S, Chau A, Matloff J, et al. The St. Jude Medical valve. *J Thorac Cardiovasc Surg* 1994;108:1010-20
- Kratz JM, Crawford Jr FA, Sade RM, et al. St. Jude Prosthesis for aortic and mitral valve replacement: a ten-year experience. *Ann Thorac Surg* 1993;56:462-68
- Masters RG, Pipe AL, Walley VM, et al. Comparative results with the St. Jude Medical and Medtronic Hall mechanical valves. *J Thorac Cardiovasc Surg* 1995;110:663-71
- Nakano K, Koyanagi H, Hashimoto A, et al. Twelve years' experience with the St. Jude Medical valve prosthesis. *Ann Thorac Surg* 1994;57:697-703
- Nistal JF, Hurle A, Revuelta JM, et al. Clinical experience with the Carbomedics valve: early results with a new bileaflet mechanical prosthesis. *J Thorac Cardiovasc Surg* 1996;112:59-68
- Peterseim DS, Cen YY, Cheruvu S, et al. Long-term outcome after biologic versus mechanical aortic valve replacement in 841 patients. *J Thorac Cardiovasc Surg* 1999;117:890-97
- Remadi JP, Bizouarn P, Baron O, et al. Mitral valve replacement with the St. Jude Medical prosthesis: a 15-year follow-up. *Ann Thorac Surg* 1998;66:762-67
- Rodler SM, Moritz A, Schreiner W, et al. Five-year follow-up after heart valve replacement with the CarboMedics bileaflet prosthesis. *Ann Thorac Surg* 1997;63:1018-25
- Skoularigis J, Essop MR, Skudicky D, et al. Frequency and severity of intravascular hemolysis after left-sided cardiac valve replacement with Medtronic Hall and St. Jude Medical prostheses and influence of prosthetic type, position, size and number. *Am J Cardiol* 1993;71:587-91
- Smith JA, Westlake GW, Mullerworth MH, et al. Excellent long-term results of cardiac valve replacement with the St. Jude Medical Valve prosthesis. *Circulation* 1993;88:49-54
- Zellner JL, Kratz JM, Crumbley III AJ, et al. Long-term experience with the St. Jude Medical valve prosthesis. *Ann Thorac Surg* 1999;68:1210-18
- David TE, et al., Correspondence: Mechanical Valves, To the Editor, *Ann Thorac Surg* 1996;62:1565-1574. See pp1567-1570.
- On-X PMA application to the U.S. Food and Drug Administration, August 2000.

Hemodynamic Comparison of Stented Pericardial, Stentless Porcine and On-X Valves

Table II
Valve Effective Orifice Area (cm²)

Size	Mitroflow	Perimount	Toronto	Freestyle	On-X
19.0	1.14 ± 0.16	0.90 ± 0.13		1.26 ± 0.27	1.5 ± 0.3
19.0		1.1 ± 0.23			
21.0				1.52 ± 0.54	1.9 ± 0.5
22.0			1.3 ± 0.7		
23.0			1.5 ± 0.5	1.77 ± 0.59	2.47 ± 0.67
25.0			1.7 ± 0.4	2.08 ± 0.62	2.7 ± 0.7
27.0			2.0 ± 0.4	2.54 ± 0.74	2.9 ± 0.7
29.0			2.4 ± 0.6		2.9 ± 0.7

Table III
Mean Gradient (mmHg)

Size	Mitroflow	Perimount	Toronto	Freestyle	On-X
19.0	14.9 ± 3.5	19.1 ± 4.9		12.1 ± 4.9	8.9 ± 3.1
19.0		18.0 ± 6.9			
21.0				9.6 ± 7.3	7.6 ± 2.96
22.0			7.3 ± 4.4		
23.0			7.4 ± 4.5	8.7 ± 7.8	6.6 ± 3.1
25.0			6.1 ± 3.1	5.9 ± 4.4	4.3 ± 2.4
27.0			4.9 ± 2.4	4.2 ± 3.0	5.6 ± 3.1
29.0			4.0 ± 2.1		5.6 ± 3.1

Table IV
Peak Gradient (mmHg)

Size	Mitroflow	Perimount	Toronto	Freestyle	On-X
19.0	20.3 ± 5.8	25.4 ± 6.4			17.1 ± 5.3
19.0		31.0 ± 11.2			
21.0					14.2 ± 5.4
22.0			18.4 ± 11.8		
23.0			15.1 ± 8.8		12.4 ± 6.2
25.0			11.6 ± 6.6		8.9 ± 4.6
27.0			9.6 ± 5.0		9.8 ± 5.3
29.0			7.2 ± 4.1		9.8 ± 5.3

** All data late post operation.*

Hemodynamic Comparison Data Sources:


1. Bach DS, David T, Yacoub M, et.al. Hemodynamics and Left Ventricular Mass Regression Following Implantation of the Toronto SPV Stentless Porcine Valve. *Amer J Cardiol* 82,1214-1219(1998)
2. Wang Z, Grainger N, Chambers J. Doppler Echocardiography in Normally Functioning Replacement Heart Valves: A Literature Review. *J Heart Valve Dis* 4,591-614(1995)
3. Dumesnil J, LeBlanc MH, Cartier P, Metraw J, Desaulniers D, Doyle D, Lemieux M, Raymond G. Hemodynamic Features of the Freestyle Aortic Bioprosthesis Compared with Stented Bioprosthesis. *Ann Thorac Surg* 66,S130-133(1998)
4. Gonzalez-Juanatey J, Garcia-Benochea J, Garcia-Acuna J, Fernandez M, Cendon A, Roman A, Pena M. The Influence of the Design on Medium to Long Term Hemodynamic Behavior of the 19 mm Pericardial Aortic Valve Prostheses. *J. Heart Valve Disease* 5;SIII:S317-S323,1996.
5. Freestyle Aortic Root Bioprosthesis: Product Performance Report, Data Current to October, 1996, Medtronic (1997)
6. On-X PMA application to the U.S. Food and Drug Administration, August 2000.

Acknowledgements

The Medical Carbon Research Institute would like to express its appreciation for the efforts of the primary clinical investigators and to its engineering staff and consultants for submitting the data reports and publications used to produce this edition of “The

On-X[®] Experience” and looks forward to their continued support.

Your comments, questions and ideas regarding “The On-X[®] Experience” should be directed to Robert More (512) 339-8000 or rmore@mcritx.com.



Headquarters:

8200 Cameron Road, Suite A-196
Austin, Texas 78754 U.S.A.
Telephone: (512) 339-8000
Facsimile: (512) 339-3636
www.mcric.com
onx@mcritx.com

Manufacturing Facilities:

8200 Cameron Road, Suite A-196
Austin, Texas U.S.A.
Telephone: (512) 339-8000
Facsimile: (512) 339-3636

Ostpassage 11
30853 Langenhagen, Germany
Telephone: 49-511-724-2820
Facsimile: 49-511-724-2817