

# CryoFreezer

## Allograft Storage Option for Hospitals



**CryoLife**<sup>®</sup>  
Life Restoring Technologies<sup>®</sup>

# Benefits to the Hospital



## Reduce Costs:

- CryoLife offers the ability to have a CryoFreezer placed in your account at **no** charge
- Avoid the return and restocking fees by storing allografts in the hospital's CryoFreezer
- Avoid the overnight shipping charges for allografts that are ordered to arrive just in time for surgery
- Save an average of \$550 plus shipping fees per return

# Benefits to the Patient

## Prepare for Unexpected Situations:

- Emergent limb salvage or infection cases at unexpected times
- The surgeon is not able to find any suitable autologous vein for bypass
- The surgeon discovers during a case that the allograft size that is needed for the patient is different from the size originally planned

## Solution:

Cryopreserved allografts are on site and available to be immediately selected from the hospital's CryoFreezer.





# CryoFreezer Options

## Small Freezer

- Storage Capacity: 24 allografts
- 3 racks for liquid storage
- Capable of dual phase storage in/ for vapor and liquid
- CryoGloves provided
- The CryoFreezer is valued at approximately \$11,360 USD (if purchased)



## Large Freezer

- Storage Capacity: 40 allografts
- 5 racks for liquid storage
- Capable of dual phase storage in/ for vapor and liquid
- CryoGloves provided
- The CryoFreezer is valued at approximately \$14,468 USD (if purchased)



## Requesting a CryoFreezer

1. Call your local CryoLife Representative or CryoLife at 1-888-427-9654
2. Your CryoLife Representative will assist you in completing the necessary forms and agreements
3. After the completion of the approval process, CryoLife will ship your hospital a freezer
4. The hospital will need to order in advance two LS-160 liquid nitrogen supply tanks to arrive a couple of days before the freezer
5. A CryoLife Representative will install the new freezer, program the CryoFreezer monitor, and train hospital personnel on the equipment

## Maintenance and Service

### Hospital Responsibilities:

- Supply the liquid nitrogen
- Check the temperature & liquid nitrogen levels of the CryoFreezer daily on the digital display and manually once a week with a measuring device

### CryoLife Responsibilities:

- Annual preventative maintenance performed by CryoLife Representative
- CryoLife Technical Support Services Representative is available 24 hours, 365 days a year at 1-800-438-8285

**For more information  
or to place an order call  
1.888.427.9654**

1655 Roberts Boulevard, NW • Kennesaw, Georgia 30144 • USA  
Tel: 770.419.3355 • 800.438.8285 • Fax: 770.590.3753  
[www.cryolife.com](http://www.cryolife.com)

All registered trademarks owned by CryoLife, Inc. or its subsidiaries.  
© 2018 CryoLife, Inc. All rights reserved.

MLENG0412.002 (02/2018)

